HUMAN CARE

Your Life. Your Way.

Altair, HeliQ, HeliQ with Power Traverse, Roomer S

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English:

Altair:

| | |
|-------------------------------------|----|
| Exploded view - Altair | 4 |
| Spare part kits | 5 |
| Opening the lift - Altair | 8 |
| Closing the lift - Altair | 9 |
| Reseting the service alarm – Altair | 9 |
| Changing the belt - Altair | 10 |
| Changing the battery - Altair | 12 |

HeliQ:

| | |
|------------------------------------|----|
| Changing the hand control - Altair | 13 |
| Exploded view - HeliQ | 14 |
| Spare part kits - HeliQ | 15 |
| Opening the lift - HeliQ | 18 |
| Closing the lift - HeliQ | 19 |
| Reseting service alarm – HeliQ | 19 |
| Changing the belt - HeliQ | 20 |
| Changing the battery - HeliQ | 21 |

Roomer S:

| | |
|--------------------------------------|----|
| Changing the hand control - HeliQ | 22 |
| Exploded view - Roomer S | 24 |
| Spare part kits - Roomer S | 25 |
| Opening the lift - Roomer S | 28 |
| Closing the lift - Roomer S | 29 |
| Reseting service alarm – Roomer S | 29 |
| Changing the hand control - Roomer S | 30 |
| Changing the belt - Roomer S | 31 |
| Changing the battery - Roomer S | 33 |

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
Exploded view - Altair

ENG

All spare parts for the overhead lifts are sold as kits. circled and pointed parts shows the kits that are available. Detailed list with included parts is shown at the next page.

| Kit number in exploded view | Part no. | Spare parts kit | Included in kit | Quantity of pieces in kit |
|-----------------------------|----------|---------------------------------------|--|---------------------------|
| 1a | 71461 | Top cover Altair - kit | Screw PT WN1452 K30x40 FZB | 2 |
| | | | Top cover Altair | 1 |
| 1b | 92937 | Bottom Cover Altair | Cover 54027 | 1 |
| | | | Label at customer specification (picture of original label required) | 1 |
| 2 | 71445 | Foot portable lift - kit | Foot portable lift | 20 |
| 3 | 71463 | Upper limit switch Altair HeliQ - kit | Belt guide Altair HeliQ | 1 |
| | | | Belt stop switch | 2 |
| | | | Belt idler shaft | 1 |
| | | | Screw MC6S M6x16 A2 | 2 |
| 4 | 71449 | Battery overhead lift - kit | CT2.3-12_Battery | 2 |
| | | | 3M VHB 5930 Double sided adhesive tape | 2 |
| 4 | 71450 | Battery overhead lift - kit | CT2.3-12_Battery | 10 |
| | | | 3M VHB 5930 Double sided adhesive tape | 10 |
| 5 | 71466 | Beltstop sensor Altair HeliQ - kit | Belt stop sensor with 350 cable | 2 |
| | | | Belt tension guard sensor 1 | 1 |
| | | | Spring Belt Tension Guard | 1 |
| | | | Washer Plastic BTG | 1 |
| | | | Washer-LB-A-4816 | 1 |
| | | | Belt tension guard shaft | 1 |
| | | | Screw MC6S M3x8 A2-70 | 3 |
| | | | Screw MC6S M6x16 A2 | 2 |
| 6 | 71474 | Circuit board Altair 150kg - kit | Printed circuit board main Altair 150kg | 1 |
| | | | Screw PT WN1451 K30x8 FZB | 2 |
| | | | Washer Plastic 4.83x3.05x2.54 | 2 |
| 6 | 71475 | Circuit board Altair 220kg - kit | Printed circuit board main Altair 220kg | 1 |
| | | | Screw PT WN1451 K30x8 FZB | 2 |
| | | | Washer Plastic 4.83x3.05x2.54 | 2 |
| 6 | 71476 | Circuit board Altair 300kg - kit | Printed circuit board main Altair 300kg | 1 |
| | | | Screw PT WN1451 K30x8 FZB | 2 |
| | | | Washer Plastic 4.83x3.05x2.54 | 2 |
| not in view | 54437 | Coin battery CR-1025_BN | Battery for main circuit board | 1 |
| 7 | 71489 | Emergency stop portable lift - kit | Emergency Stop 16 | 1 |

| Kit number in exploded view | Part no. | Spare parts kit | Included in kit | Quantity of pieces in kit |
|-----------------------------|----------|--|--|---------------------------|
| 8a | 71483 | Lift belt ceiling lift SWL 150kg/330lbs and 220kg/440lbs | Lift belt for Altair, Roomer S and HeliQ with lifting length 2120 mm | 1 |
| 8b | 71494 | Lift belt ceiling lift SWL 300kg/660lbs - kit | Lift belt for ceiling lift Altair, Roomer S, HeliQ with lifting length 1720 mm | 1 |
| 9 | 71467 | Belt roll Altair HeliQ - kit | Screw MC6S M6x16 A2 | 6 |
| | | | Belt roll Altair HeliQ | 1 |
| | | | Shaft belt roll Altair HeliQ | 1 |
| | | | Washer roll | 1 |
| 10a | 71455 | Spring telescopic hanger bar ceiling lift - kit | Spring | 10 |
| 10b | 72803 | Locking clips for hanger bar | Clips for hanger bars 72754, 72755, 72757, 55751, 55764, 92938 | 10 |
| 10c | 92764 | Shaft for hanger bar | For assembly of hanger bars 755751, 55764 | 1 |
| 11a | 71456 | Spring stainless steel leg hook Altair, Roomer S - kit | Spring | 5 |
| 11b | 92974 | Spring coated leg hook - 5 pack | Spring for leg hook | 5 |
| 12 | 71457 | Hand control Altair Roomer S - kit | Hand control Altair Roomer S without covers | 1 |
| | | | Strain Relief Clamp cc12 | 1 |
| | | | Screw PT WN1452 K30x10 FZB | 2 |
| 13 | 71462 | Screw holder Altair HeliQ - kit | Screw Holder | 10 |
| 14 | 71813 | Secondary control panel Roomer S Altair HeliQ | Membrane panel Roomer S Altair HeliQ | 1 |
| not in view | 71470 | Charger Europe, overhead lift - kit | Battery Charger | 1 |
| | | | Charger Cable Europe | 1 |
| | | | Hand control holder | 1 |
| | | | Ring strap for holder | 1 |
| not in view | 71471 | Charger Great Britain, overhead lift - kit | Battery Charger | 1 |
| | | | Cable Charger Great Britain | 1 |
| | | | Hand control holder | 1 |
| | | | Ring strap for holder | 1 |
| not in view | 71472 | Charger Japan, overhead lift - kit | Battery Charger | 1 |
| | | | Cable Charger Japan | 1 |
| | | | Hand control holder | 1 |
| | | | Ring strap for holder | 1 |
| not in view | 71473 | Charger North America, overhead lift - kit | Battery Charger | 1 |
| | | | Cable Charger North America | 1 |
| | | | Hand control holder | 1 |

| Kit number in exploded view | Part no. | Spare parts kit | Included in kit | Quantity of pieces in kit |
|------------------------------------|-----------------|---|---------------------------------------|----------------------------------|
| | | | Ring strap for holder | 1 |
| not in view | 71491 | Charger Australia overhead lift - kit | Battery charger | 1 |
| | | | Cable charger Australia | 1 |
| | | | Hand control holder | 1 |
| | | | Ring strap for holder | 1 |
| not in view | 90407 | Improvement color RAL9016 | | 1 |
| not in view | 71811 | Cover screw kit | Cover screw | 10 |
| not in view | 92972 | Secondary control panel RoomerS, Altair HeliQ | Membrane panel RoomerS, Altair, HeliQ | 1 |

Opening the lift - Altair

ENG


1. Activate the emergency stop.


2. Turn the lift upside down. Be careful with the belt stop sensors.


3. Remove the rubber feet closest to the leg hook on the Altair lift.


4. Untighten the two screws with a torx screwdriver (TX10) holding the covers together (about 6-7 laps or until it's loosened). The screws are located at the bottom of the holes for the rubber feet/plastic.

5. Turn the lift around and remove the cover by angle it a bit towards the hanger bar until you hear a "click".

6. Always disconnect the battery cable when the cover is opened.

1. Make sure that all cables are connected as they were when opening the lift.
2. Make sure no cables are pinched between the covers.
3. Follow the steps 1-6 in the “opening the lift chapter” in reverse.

Don't forget the rubber feet!

Remember to attach all cables before closing the lift!

Reseting the service alarm – Altair

After the service is finished the service alarm should be reset.

- Connect the lift to a computer thru “Data and service application”. In the software, click service and log in with your login details. Click on the button “Reset service alarm”.


The screenshot shows the 'HUMAN CARE Data- och serviceapplikation' interface. The 'SERVICE' tab is active. A button labeled 'ÅTERSTÄLL SERVICEALARM' is circled in red. Below it, there is a table showing service history with columns for Datum, Tid, and Beskrivning.

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
The belt can be changed without removing the belt outlet.

Follow the steps 1-5 in the “opening the lift chapter”.


1. Re-activate the emergency stop.
2. Unscrew the Opto-switch for the belt tension guard and turn it 90° in order to be able to freely operate the lift without weight.


3. Drag the belt out all the way from the belt roll.

4. Create a big loop.

5. Use for example a spatula to create a loop on the other side.


6. Press “up” button on the hand control to turn the belt roll and grab the loop.

7. Use the spatula again to get the end of the belt with the small pin out.


8. Return the belt roll with the up button on the hand control. Grab the belt.

9. Remove the pin, the belt is now loose.

10. To insert the new belt, do the steps in reverse.


After the new belt is assembled. Make sure to position the opto-switch to it's starting position!


The battery pack consists of two 12 V, 2.3 Ah, lead type batteries. The lifetime is up to 5 years depending on charging and usage.

Replacing the battery:


1. Follow the steps 1-6 in the “opening the lift chapter”.
2. Try to not damage the warning label on top of the battery, for example by folding up the ends of the label.


3. Unscrew the two metal brackets from the chassis by using a 2,5 mm Allen key.

4. Lift up the battery pack from the lift.


5. When the batteries is removed from the lift, disconnect all four cables.


6. Fit the new batteries and re-connect the cables. Make sure the batteries are connected the same way as before the change.
7. Secure the battery cables between the battery (outside the steel brackets) and the lift cover.
8. Follow the steps 1-3 in reverse to complete the battery change.


1. Follow the steps 1-6 in the “opening the lift chapter”.
2. Remove the connector on the hand control cable from the circuit board (PCB) by pressing on the small plastic lock on the left side.



4. Remove the rubber seal from the cover.

5. Pull the cable and the connector out through the hole.

3. Remove the two screws from the cable holder with a torx screwdriver 10 mm (T10).


6. Repeat the steps in reverse in order to connect the new hand control.

If the hand control is damaged/broken, the cable can be cut to make it easier to remove the cable.

Exploded view - HeliQ

ENG

All spare parts for the overhead lifts are sold as a kits. circled and pointed parts shows the kits that are available. Detailed list with included parts is shown at the next page.

| Kit number in exploded view | Part no. | Spare parts kit | Included in kit | Quantity of pieces in kit |
|-----------------------------|----------|---------------------------------------|---|---------------------------|
| 1a | 71464 | Bottom cover HeliQ - kit | Screw PT WN1452 K30x40 FZB | 2 |
| | | | Bottom cover | 1 |
| 1b | 92971 | Top cover HeliQ | Cover 54028 | 1 |
| | | | Label at customer specifications (picture of original label required) | 1 |
| 2 | 71465 | Cover plug HeliQ - kit | Cover Plug 9.3 | 20 |
| 3 | 71463 | Upper limit switch Altair HeliQ - kit | Belt guide Altair HeliQ | 1 |
| | | | Belt stop switch | 2 |
| | | | Belt idler shaft | 1 |
| | | | Screw MC6S M6x16 A2 | 2 |
| 4 | 71449 | Battery overhead lift - kit | CT2.3-12_Battery | 2 |
| | | | 3M VHB 5930 Double sided adhesive tape | 2 |
| 4 | 71450 | Battery overhead lift - kit | CT2.3-12_Battery | 10 |
| | | | 3M VHB 5930 Double sided adhesive tape | 10 |
| 5 | 71466 | Beltstop sensor Altair HeliQ - kit | Belt stop sensor with 350 cable | 2 |
| | | | Belt tension guard sensor 1 | 1 |
| | | | Spring Belt Tension Guard | 1 |
| | | | Washer Plastic BTG | 1 |
| | | | Washer-LB-A-4816 | 1 |
| | | | Belt tension guard shaft | 1 |
| | | | Screw MC6S M3x8 A2-70 | 3 |
| | | | Screw MC6S M6x16 A2 | 2 |
| | | | Belt glider 1 pc. | 1 |
| 6 | 71477 | Circuit board HeliQ 150kg - kit | Printed circuit board main HeliQ 150kg | 1 |
| | | | Screw PT WN1451 K30x8 FZB | 2 |
| | | | Washer Plastic 4.83x3.05x2.54 | 2 |
| 6 | 71478 | Circuit board HeliQ 220kg - kit | Printed circuit board main HeliQ 220kg | 1 |
| | | | Screw PT WN1451 K30x8 FZB | 2 |
| | | | Washer Plastic 4.83x3.05x2.54 | 2 |
| 6 | 71479 | Circuit board HeliQ 300kg - kit | Printed circuit board main HeliQ 300kg | 1 |
| | | | Screw PT WN1451 K30x8 FZB | 2 |
| | | | Washer Plastic 4.83x3.05x2.54 | 2 |
| not in view | 54437 | Coin battery CR-1025_BN | Battery coin for main circuit board | 1 |

| Kit number in exploded view | Part no. | Spare parts kit | Included in kit | Quantity of pieces in kit |
|-----------------------------|----------|--|--|---------------------------|
| 7 | 71490 | Emergency stop HeliQ - kit | Emergency Stop HeliQ | 1 |
| 8a | 71483 | Lift belt ceiling lift SWL 150kg/330lbs and 220kg/440lbs | Lift belt for Altair, Roomer S and HeliQ with lifting length 2120mm | 1 |
| 8b | 71494 | Lift belt ceiling lift SWL 300kg/660lbs - kit | Lift belt for ceiling lift Altair, Roomer S, HeliQ with SWL, lifting length 1720mm | 1 |
| 9 | 71467 | Belt roll Altair HeliQ - kit | Screw MC6S M6x16 A2 | 6 |
| | | | Belt Roll Altair HeliQ | 1 |
| | | | Shaft Belt Roll Altair HeliQ | 1 |
| | | | Washer roll | 1 |
| 10a | 71455 | Spring telescopic hanger bar ceiling lift - kit | Spring | 10 |
| 10b | 72803 | Locking clips for hanger bar | Clips for hanger bars 72754, 72755, 72757, 55751, 55764, 92938 | 10 |
| 11 | 71458 | Hand control HeliQ - kit | Hand control 4-button HeliQ without covers | 1 |
| | | | Strain Relief Clamp cc12 | 1 |
| | | | Screw PT WN1452 K30x10 FZB | 2 |
| 12 | 71459 | Hand control HeliQ Power traverse - kit | Hand control 7-button HeliQ power traverse | 1 |
| | | | Strain Relief Clamp cc12 | 1 |
| | | | Screw PT WN1452 K30x10 FZB | 2 |
| 13 | 71462 | Screw holder Altair HeliQ - kit | Screw Holder | 10 |
| 14 | 71468 | Top plate screw HeliQ - kit | Screw MC6S M6x16 A2-70 | 20 |
| 15 | 71469 | Pin and locking pin HeliQ - kit | Locking needle | 2 |
| | | | Split pin CBH 20x25 fzb | 2 |
| not in view | 71470 | Charger Europe, overhead lift - kit | Battery Charger | 1 |
| | | | Cable Charger Europe | 1 |
| | | | Hand control holder | 1 |
| | | | Ring strap for holder | 1 |
| not in view | 71471 | Charger Great Britain, overhead lift - kit | Battery Charger | 1 |
| | | | Cable Charger Great Britain | 1 |
| | | | Hand control holder | 1 |
| | | | Ring strap for holder | 1 |
| not in view | 71472 | Charger Japan, overhead lift - kit | Battery Charger | 1 |
| | | | Cable Charger Japan | 1 |
| | | | Hand control holder | 1 |
| | | | Ring strap for holder | 1 |
| not in view | 71473 | Charger North America, overhead lift - kit | Battery Charger | 1 |
| | | | Cable Charger North America | 1 |
| | | | Hand control holder | 1 |

| Kit number in exploded view | Part no. | Spare parts kit | Included in kit | Quantity of pieces in kit |
|-----------------------------|----------|--|--------------------------------------|---------------------------|
| | | | Ring strap for holder | 1 |
| not in view | 71491 | Charger Australia, overhead lift - kit | Battery Charger | 1 |
| | | | Cable Charger Australia | 1 |
| | | | Hand control holder | 1 |
| | | | Ring strap for holder | 1 |
| not in view | 90407 | Improvement color RAL9016 | | 1 |
| not in view | 71811 | Cover screw kit | Cover screw | 10 |
| not in view | 92972 | Secondary control panel RoomerS, Altair, HeliQ | Membrane panel Roomer S Altair HeliQ | 1 |

Opening the lift - HeliQ

ENG


1. Activate the emergency stop.

2. Remove the top plate screws using a 5 mm Allen key. Remove the top plate.

3. Remove the two cover plugs.


4. Untighten the two screws with a torx screwdriver (TX10) holding the covers together.

5. Turn the lift around and remove the cover by angle it a bit towards the opposite side of the control panel until you hear a "click".

6. Always disconnect the battery cable when the cover is opened.

1. Make sure that all cables are connected as they were when opening the lift.
2. Make sure no cables are pinched between the covers.
3. Follow the steps 1-6 in the “opening the lift chapter” in reverse.

Don't forget the cover plugs!

Remember to attach all cables before closing the lift!

Reseting service alarm – HeliQ

After the service is finished the service alarm should be reset.


- Connect the lift to a computer thru “Data and service application”. In the software, click service and log in with your login details. Click on the button “Reset service alarm”.

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
The belt can be changed without removing the belt outlet.

Follow the steps 1-5 in the “opening the lift chapter”.


1. Re-activate the emergency stop.
2. Unscrew the Opto-switch for the belt tension guard and turn it 90° in order to be able to freely operate the lift without weight.


3. Drag the belt out all the way from the belt stop.

4. Create a big loop.

5. Use for example a spatula to create a loop on the other side.


6. Press “up” button on the hand control to turn the belt roll and grab the loop.

7. Use the spatula again to get the end of the belt with the small pin out.


8. Return the belt roll with help of the hand control (up-button). Grab the belt.

9. Remove the pin and the belt is loose.

9. To insert the new belt, do the steps in reverse.


After the new belt is assembled. Make sure to position the opto-switch to it's starting position!


The battery pack consists of two 12 V, 2.3 Ah, lead type batteries. The lifetime is up to 5 years depending on charging and usage.

Replacing the battery:


1. Follow the steps 1-6 in the “opening the lift chapter”.
2. Try to not damage the warning label on top of the battery, for example by folding up the ends of the label.


3. Unscrew the two metal brackets from the chassis by using a 2,5 mm Allen key.

4. Lift up the battery pack from the lift.


4. When the batteries is removed from the lift, disconnect all four cables.



5. Fit the new batteries and re-connect the cables. Make sure the batteries are connected the same way as before the change.
6. Secure the battery cables between the battery (outside the steel brackets) and the lift cover.
7. Follow the steps 1-3 in reverse to complete the battery change.


1. Follow the steps 1-6 in the “opening the lift chapter”.
2. Remove the connector on the hand control cable from the circuit board (PCB) by pressing on the small plastic lock on the left side.


3. Remove the two screws from the cable holder with a torx screwdriver 10 mm (T10).

4. Remove the rubber seal from the cover.

5. Pull the cable and the connector out through the hole.


6. Repeat the steps in reverse in order to connect the new hand control.

If the hand control is damaged/broken the cable can be cut for easier removing of the hand control and cable.

Exploded view - Roomer S

ENG

All spare parts for the overhead lifts are sold as a kits. circled and pointed parts shows the kits that are available. Detailed list with included parts is shown at the next page.

| Kit number (exploded view) | Part no. | Spare parts kit | Included in kit | Quantity of pieces in kit |
|----------------------------|----------|-----------------------------------|---|---------------------------|
| 1a | 71460 | Top cover Roomer S - kit | Screw PT WN1452 K30x40 fzb | 2 |
| | | | Top cover Roomer S | 1 |
| 1b | 92936 | Bottom Cover Roomer S kit | Cover 54055 | 1 |
| | | | Label at customer specifications (picture of original label required) | 1 |
| 2 | 71445 | Foot Portable lift - kit | Foot portable lift | 20 |
| 3 | 71446 | Slider Roomer S - kit | Screw MC6LS M6x16 | 3 |
| | | | Screw MC6S M6x16 A2 | 4 |
| | | | Slider with sensor | 2 |
| 4 | 71447 | Switch Roomer S - kit | Roomer S switch lever | 2 |
| | | | Shaft Slider | 1 |
| | | | Screw MC6S M3x12 A2 | 4 |
| | | | 10015-N08 Mekanex | 2 |
| | | | Slide bearing 8x9x3 | 2 |
| 5 | 71448 | Upper limit switch Roomer S - kit | Belt guide Roomer S | 1 |
| | | | Belt stop switch Roomer S | 2 |
| | | | Belt idler shaft Roomer S | 2 |
| | | | Screw MC6S M6x16 A2 | 4 |
| 6 | 71449 | Battery overhead lift - kit | CT2.3-12_Battery | 2 |
| | | | 3M VHB 5930 Double sided adhesive tape | 2 |
| 6 | 71450 | Battery overhead lift - kit | CT2.3-12_Battery | 10 |
| | | | 3M VHB 5930 Double sided adhesive tape | 10 |
| 7a | 71451 | Beltstop sensor Roomer S - kit | Belt stop sensor with 350 cable | 1 |
| | | | Belt stop sensor with 400 cable | 1 |
| | | | Belt tension guard sensor 2 | 1 |
| | | | Belt tension guard sensor 1 | 1 |
| | | | Spring belt tension guard | 4 |
| | | | Washer Plastic BTG | 4 |
| | | | Washer-LB-A-4816 | 2 |
| | | | Belt tension guard shaft Belt tension guard shaft | 2 |
| | | | Screw MC6S M3x8 A2-70 | 4 |
| | | | Screw MC6S M6x16 A2 | 4 |
| | | | Belt glider | 2 |
| | | | Cable ties 2,5x100 | 3 |
| | | | 7b | 71551 |

| Kit number in exploded view | Part no. | Spare parts kit | Included in kit | Quantity of pieces in kit |
|-----------------------------|----------|--|--|---------------------------|
| 8 | 71452 | Circuit board Roomer S - kit | Printed circuit board main Roomer S 220 kg | 1 |
| | | | Screw MC6S M3x8 A2-70 | 3 |
| not in view | 54437 | Coin battery CR-1025_BN | Battery coin for main circuit board | 1 |
| 9 | 71489 | Emergency stop portable lifts - kit | Emergency Stop 16 | 1 |
| 10a | 71483 | Lift belt ceiling lift SWL 150kg/330lbs and 220kg/440lbs | Lift belt for Altair, Roomer S and HeliQ with lifting length 2120mm | 1 |
| 10b | 71494 | Lift belt ceiling lift SWL 300kg/660lbs - kit | Lift belt for ceiling lift Altair, Roomer S, HeliQ with SWL, lifting length 1720mm | 1 |
| 11a | 71455 | Spring telescopic hanger bar ceiling lift - kit | Spring | 10 |
| 11b | 72803 | Locking clips for hanger bar | Clips for hanger bars 72754, 72755, 72757, 55751, 55764, 92938 | 10 |
| 11c | 92764 | Shaft hanger bar Altair, Roomer S | For assembly of hanger bar 55751 and 55764 | 1 |
| 12 | 71453 | Belt roll clock wise Roomer S -kit | Screw MC6S M6x16 A2 | 4 |
| | | | Belt roll clock wise | 1 |
| | | | Bushing for belt roll shaft | 2 |
| | | | Shaft Belt Roll | 1 |
| 13 | 71454 | Belt roll counter clock wise Roomer S - kit | Screw MC6S M6x16 A2 | 4 |
| | | | Belt roll counter clock wise | 1 |
| | | | Bushing for belt roll shaft | 2 |
| | | | Shaft Belt Roll | 1 |
| 14a | 71456 | Spring stainless steel leg hook Altair Roomer S - kit | Spring | 5 |
| 14b | 92974 | Spring coated leg hook - 5 pack | Spring for leg hook | 5 |
| 15 | 71457 | Hand control Altair Roomer S - kit | Hand control portable lifts without cover | 1 |
| | | | Strain Relief Clamp cc12 | 1 |
| | | | Screw PT WN1452 K30x10 FZB | 2 |
| not in view | 71470 | Charger Europe, overhead lift - kit | Battery Charger | 1 |
| | | | Cable Charger Europe | 1 |
| | | | Hand control holder | 1 |
| | | | Ring strap for holder | 1 |
| not in view | 71471 | Charger Great Britain, overhead lift - kit | Battery Charger | 1 |
| | | | Cable Charger Great Britain | 1 |
| | | | Hand control holder | 1 |

| Kit number in exploded view | Part no. | Spare parts kit | Included in kit | Quantity of pieces in kit |
|-----------------------------|----------|--|--------------------------------------|---------------------------|
| | | | Ring strap for holder | 1 |
| not in view | 71472 | Charger Japan, overhead lift - kit | Battery Charger | 1 |
| | | | Cable Charger Japan | 1 |
| | | | Hand control holder | 1 |
| | | | Ring strap for holder | 1 |
| not in view | 71473 | Charger North America, overhead lift - kit | Battery Charger | 1 |
| | | | Cable Charger North America | 1 |
| | | | Hand control holder | 1 |
| | | | Ring strap for holder | 1 |
| not in view | 71491 | Charger Australia, overhead lift - kit | Battery Charger | 1 |
| | | | Cable Charger Australia | 1 |
| | | | Hand control holder | 1 |
| | | | Ring strap for holder | 1 |
| not in view | 90407 | Improvement color RAL9016 | | 1 |
| not in view | 71811 | Cover screw kit | Cover screw | 10 |
| not in view | 92972 | Secondary control panel RoomerS, Altair, HeliQ | Membrane panel Roomer S Altair HeliQ | |
| not in view | 92973 | Battery holder kit | 54905 Battery holder RoomerS | 1 |
| | | | 55629 Heat shrink tubing | 1 |

Opening the lift - Roomer S

ENG

1. Activate the emergency stop.


2. Turn the lift upside down. Be careful with the belt stop sensor.

3. Remove the rubber feet closest to the hanger bar on Roomer S lift.


4. Untighten the two screws (TX10) holding the covers together. The screws are located at the bottom of the holes for the rubber feet.

5. Turn the lift around and remove the cover by angle it a bit towards the opposite side of the control panel until you hear a "click".

6. Always disconnect the batteries when the cover is opened.

1. Make sure that all cables are connected as they were before.
2. Make sure no cables are pinched.
3. Follow the steps 1-6 in the “opening the lift chapter” in reverse.

Remember to attach the battery cable before closing the lift!

Don't forget the rubber feet!

Reseting service alarm – Roomer S

After the service is finished the service alarm should be reset.

- Connect the lift to a computer thru “Data and service application”. In the software, click service and log in with your login details. Click on the button “Reset service alarm”.



The screenshot shows the 'Data- och serviceapplikation' interface for HUMAN CARE. The top navigation bar includes 'START', 'ÖVERSIKT', 'STATISTIK', and 'SERVICE'. The 'SERVICE' tab is active. Below the navigation bar, there are three main sections: 'ÅTERSTÄLL SERVICEALARM' (highlighted with a red circle), 'UPPDATERA TIDPUNKT', and 'STÄLL IN REALTIDSKLOCKA'. The 'SERVICE' section displays 'Senaste service' (2015-04-08) and 'Nästa service tillfälle' (2016-04-07). The 'FELRAPPORTER' section shows a table of error reports.

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
Changing the hand control - Roomer S

ENG

1. Follow the steps 1-6 in the “opening the lift chapter”.
2. Remove the connector on the hand control cable from the circuit board (PCB) by pressing on the small plastic lock on the top.


3. Remove the rubber seal from the cover.


4. Use for example a screwdriver to flip over

the cable to the left side like in pictures below.

5. Pull the cable and the connector out through the hole. If the hand control is damaged, you can cut the cable and it will be easier to pull out.


6. Repeat the steps in reverse in order to connect the new hand control.

Changing the belt - Rooter S


ENG

The belt can be changed without removing the belt outlet.


1. Drag the unlocked belt out all the way from the belt roll counter.


2. Create a quite big loop.

3. Use for example a spatula to create a loop on the other side.

5. Use the roomer S switch or the hand control to turn the belt roll and grab the loop.

6. Use the spatula again to get the end of the belt with the small pin out.


7. Return the belt roll with help of the Rooter S switch or the hand control. Grab the belt.

Note! Do not turn the Rooter Switch without holding the belt roll since the spring will unwind. If you drop the belt roll, rewind the belt roll 9,5 turns from the belt roll slot with a metal pin.


7. Remove the pin and the belt is loose.

8. If you need to remove both belts, do the same on the opposite side.

9. Insert the new belt/belts through the belt outlet as shown in the picture below.

Note! Make sure the black tags on both belts are facing each other.

Note! Make sure the belts are on the correct side of the Belt tension guard.


- The belt should enter the belt roll through the slot with the metal pin.
- The end of the belt should be fastened with the metal pin in the slot next to the slot where it entered the belt roll.
- The pin should be resting in the slot as much as possible. Pull the belt to make the pin fit in the slot.

The battery pack consists of two 12 V, 2.3 Ah, lead type batteries. The lifetime is up to 5 years depending on charging and usage.

Replacing the battery:


1. Follow the steps 1-6 in the “opening the lift chapter”.

Note! Make sure to remove the battery cable that connects to the circuit board (not the battery cables that is placed on the battery).


2. Unscrew the metal bracket from the chassis by using a 2,5 mm Allen key.

Warning! Do not remove the battery cables show in the picture below! They should stay on until the screw is loosened to avoid short circuit.


3. Now - disconnect the battery cables located on the battery.

4. Lift up the battery pack from the lift.


5. When the batteries is removed from the lift, disconnect all four cables.

5. Fit the new batteries and re-connect the cables. Make sure the batteries are connected the same way as before. The batteries connectors from the circuit board and the battery terminals are color coded. Cables from emergency stop can be connected in any order.
6. Secure the battery cables between the battery (outside the steel brackets) and the lift cover.
7. Follow the steps 1-5 in reverse to complete the battery change.

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| Kitnummer i sprängskiss | Art. nr. | Reservdelskit | Inkluderat i kit | Antal delar i kit |
|-------------------------|----------|---|---|-------------------|
| 1a | 71460 | Övre kåpa Roomer S - kit | Screw PT WN1452 K30x40 fzb | 2 |
| | | | Top cover Roomer S | 1 |
| 1b | 92936 | Undre kåpa Roomer S kit | Kåpa 54055 | 1 |
| | | | Etikett vid kundspecifikation (bild av original etikett nödvändig) | 1 |
| 2 | 71445 | Fötter portabel lyft - kit | Foot portable lift | 20 |
| 3 | 71446 | Glidplatta Roomer S - kit | Screw MC6LS M6x16 | 3 |
| | | | Screw MC6S M6x16 A2 | 4 |
| | | | Slider with sensor | 2 |
| 4 | 71447 | Spak Roomer S - kit | Roomer S switch lever | 2 |
| | | | Shaft Slider | 1 |
| | | | Screw MC6S M3x12 A2 | 4 |
| | | | 10015-N08 Mekanex | 2 |
| | | | Slide bearing 8x9x3 | 2 |
| 5 | 71448 | Övre brytare Roomer S - kit | Belt guide Roomer S | 1 |
| | | | Belt stop switch Roomer S | 2 |
| | | | Belt idler shaft Roomer S | 2 |
| | | | Screw MC6S M6x16 A2 | 4 |
| 6 | 71449 | Batteri taklyft - kit 2 st | CT2.3-12_Battery | 2 |
| | | | 3M VHB 5930 Double sided adhesive tape | 2 |
| 6 | 71450 | Batteri taklyft - kit 10 st | CT2.3-12_Battery | 10 |
| | | | 3M VHB 5930 Double sided adhesive tape | 10 |
| 7 | 71451 | Sensorer bandstyrning Roomer S - kit | Beltstop sensor with 350 cable | 1 |
| | | | Beltstop sensor with 400 cable | 1 |
| | | | Belttension guard sensor 2 | 1 |
| | | | Belttension guard sensor 1 | 1 |
| | | | Spring belt tension guard | 4 |
| | | | Washer Plastic BTG | 4 |
| | | | Washer-LB-A-4816 | 2 |
| | | | Belt tension guard shaft | 2 |
| | | | Screw MC6S M3x8 A2-70 | 4 |
| | | | Screw MC6S M6x16 A2 | 4 |
| Cable ties 2,5x100 | 3 | | | |

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