HUMAN CARE

Your Life. Your Way.

Altair, HeliQ, HeliQ with Power Traverse, Roomer S

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English:

Altair:

Exploded view - Altair	4
Spare part kits	5
Opening the lift - Altair	8
Closing the lift - Altair	9
Reseting the service alarm – Altair	9
Changing the belt - Altair	10
Changing the battery - Altair	12

HeliQ:

Changing the hand control - Altair	13
Exploded view - HeliQ	14
Spare part kits - HeliQ	15
Opening the lift - HeliQ	18
Closing the lift - HeliQ	19
Reseting service alarm – HeliQ	19
Changing the belt - HeliQ	20
Changing the battery - HeliQ	21

Roomer S:

Changing the hand control - HeliQ	22
Exploded view - Roomer S	24
Spare part kits - Roomer S	25
Opening the lift - Roomer S	28
Closing the lift - Roomer S	29
Reseting service alarm – Roomer S	29
Changing the hand control - Roomer S	30
Changing the belt - Roomer S	31
Changing the battery - Roomer S	33

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
Exploded view - Altair

ENG

All spare parts for the overhead lifts are sold as kits. circled and pointed parts shows the kits that are available. Detailed list with included parts is shown at the next page.

Kit number in exploded view	Part no.	Spare parts kit	Included in kit	Quantity of pieces in kit
1a	71461	Top cover Altair - kit	Screw PT WN1452 K30x40 FZB	2
			Top cover Altair	1
1b	92937	Bottom Cover Altair	Cover 54027	1
			Label at customer specification (picture of original label required)	1
2	71445	Foot portable lift - kit	Foot portable lift	20
3	71463	Upper limit switch Altair HeliQ - kit	Belt guide Altair HeliQ	1
			Belt stop switch	2
			Belt idler shaft	1
			Screw MC6S M6x16 A2	2
4	94356	Battery ceiling lift - kit	Battery 12V, 2.9Ah	2
			3M VHB 5930 Double sided adhesive tape	2
5	71466	Beltstop sensor Altair HeliQ - kit	Belt stop sensor with 350 cable	2
			Belt tension guard sensor 1	1
			Spring Belt Tension Guard	1
			Washer Plastic BTG	1
			Washer-LB-A-4816	1
			Belt tension guard shaft	1
			Screw MC6S M3x8 A2-70	3
			Screw MC6S M6x16 A2	2
6	71474	Circuit board Altair 150kg - kit	Printed circuit board main Altair 150kg	1
			Screw PT WN1451 K30x8 FZB	2
			Washer Plastic 4.83x3.05x2.54	2
6	71475	Circuit board Altair 220kg - kit	Printed circuit board main Altair 220kg	1
			Screw PT WN1451 K30x8 FZB	2
			Washer Plastic 4.83x3.05x2.54	2
6	71476	Circuit board Altair 300kg - kit	Printed circuit board main Altair 300kg	1
			Screw PT WN1451 K30x8 FZB	2
			Washer Plastic 4.83x3.05x2.54	2
not in view	54437	Coin battery CR-1025_BN	Battery for main circuit board	1
7	71489	Emergency stop portable lift - kit	Emergency Stop 16	1
8a	71483	Lift belt ceiling lift SWL 150kg/330lbs and 220kg/440lbs	Lift belt for Altair, Roomer S and HeliQ with lifting length 2120 mm	1

Kit number in exploded view	Part no.	Spare parts kit	Included in kit	Quantity of pieces in kit
8b	71494	Lift belt ceiling lift SWL 300kg/660lbs - kit	Lift belt for ceiling lift Altair, Roomer S, HeliQ with lifting length 1720 mm	1
9	71467	Belt roll Altair HeliQ - kit	Screw MC6S M6x16 A2	6
			Belt roll Altair HeliQ	1
			Shaft belt roll Altair HeliQ	1
			Washer roll	1
10a	71455	Spring telescopic hanger bar ceiling lift - kit	Spring	10
10b	72803	Locking clips for hanger bar	Clips for hanger bars 72754, 72755, 72757, 55751, 55764, 92938	10
10c	92764	Shaft for hanger bar	For assembly of hanger bars 755751, 55764	1
11a	71456	Spring stainless steel leg hook Altair, Roomer S - kit	Spring	5
11b	92974	Spring coated leg hook - 5 pack	Spring for leg hook	5
12	71457	Hand control Altair Roomer S - kit	Hand control Altair Roomer S without covers	1
			Strain Relief Clamp cc12	1
			Screw PT WN1452 K30x10 FZB	2
13	71462	Screw holder Altair HeliQ - kit	Screw Holder	10
14	71813	Secondary control panel Roomer S Altair HeliQ	Membrane panel Roomer S Altair HeliQ	1
not in view	71470	Charger Europe, overhead lift - kit	Battery Charger	1
			Charger Cable Europe	1
			Hand control holder	1
			Ring strap for holder	1
not in view	71471	Charger Great Britain, overhead lift - kit	Battery Charger	1
			Cable Charger Great Britain	1
			Hand control holder	1
			Ring strap for holder	1
not in view	71472	Charger Japan, overhead lift - kit	Battery Charger	1
			Cable Charger Japan	1
			Hand control holder	1
			Ring strap for holder	1
not in view	71473	Charger North America, overhead lift - kit	Battery Charger	1
			Cable Charger North America	1
			Hand control holder	1
			Ring strap for holder	1

Kit number in exploded view	Part no.	Spare parts kit	Included in kit	Quantity of pieces in kit
not in view	71491	Charger Australia overhead lift - kit	Battery charger	1
			Cable charger Australia	1
			Hand control holder	1
			Ring strap for holder	1
not in view	90407	Improvement color RAL9016		1
not in view	71811	Cover screw kit	Cover screw	10
not in view	92972	Secondary control panel RoomerS, Altair HeliQ	Membrane panel RoomerS, Altair, HeliQ	1


Opening the lift - Altair

ENG


1. Activate the emergency stop.


2. Turn the lift upside down. Be careful with the belt stop sensors.


3. Remove the rubber feet closest to the leg hook on the Altair lift.


4. Untighten the two screws with a torx screwdriver (TX10) holding the covers together (about 6-7 laps or until it's loosened). The screws are located at the bottom of the holes for the rubber feet/plastic.

5. Turn the lift around and remove the cover by angle it a bit towards the hanger bar until you hear a "click".

6. Always disconnect the battery cable when the cover is opened.

1. Make sure that all cables are connected as they were when opening the lift.
2. Make sure no cables are pinched between the covers.
3. Follow the steps 1-6 in the “opening the lift chapter” in reverse.


Don't forget the rubber feet!

Remember to attach all cables before closing the lift!

Reseting the service alarm – Altair

After the service is finished the service alarm should be reset.

- Connect the lift to a computer thru “Data and service application”. In the software, click service and log in with your login details. Click on the button “Reset service alarm”.


The screenshot shows the 'HUMAN CARE Data- och serviceapplikation' interface. The 'SERVICE' tab is selected. A red circle highlights the 'ÅTERSTÄLL SERVICEALARM' button. Below the button, there are sections for 'SERVICER' and 'FELRAPPORTER'.

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
The belt can be changed without removing the belt outlet.

Follow the steps 1-5 in the “opening the lift chapter”.


1. Re-activate the emergency stop.
2. Unscrew the Opto-switch for the belt tension guard and turn it 90° in order to be able to freely operate the lift without weight.


3. Drag the belt out all the way from the belt roll.

4. Create a big loop.

5. Use for example a spatula to create a loop on the other side.


6. Press “up” button on the hand control to turn the belt roll and grab the loop.

7. Use the spatula again to get the end of the belt with the small pin out.


8. Return the belt roll with the up button on the hand control. Grab the belt.

9. Remove the pin, the belt is now loose.

10. To insert the new belt, do the steps in reverse.


After the new belt is assembled. Make sure to position the opto-switch to it's starting position!


The battery pack consists of two 12 V, 2.9 Ah, lead type batteries. The lifetime is up to 5 years depending on charging and usage.

Replacing the battery:


1. Follow the steps 1-6 in the “opening the lift chapter”.
2. Try to not damage the warning label on top of the battery, for example by folding up the ends of the label.


3. Unscrew the two metal brackets from the chassis by using a 2,5 mm Allen key.

4. Lift up the battery pack from the lift.


5. When the batteries is removed from the lift, disconnect all four cables.


6. Fit the new batteries and re-connect the cables. Make sure the batteries are connected the same way as before the change.
7. Secure the battery cables between the battery (outside the steel brackets) and the lift cover.
8. Follow the steps 1-3 in reverse to complete the battery change.


1. Follow the steps 1-6 in the “opening the lift chapter”.
2. Remove the connector on the hand control cable from the circuit board (PCB) by pressing on the small plastic lock on the left side.



4. Remove the rubber seal from the cover.

5. Pull the cable and the connector out through the hole.

3. Remove the two screws from the cable holder with a torx screwdriver 10 mm (T10).


6. Repeat the steps in reverse in order to connect the new hand control.

If the hand control is damaged/broken, the cable can be cut to make it easier to remove the cable.

Exploded view - HeliQ

ENG

All spare parts for the overhead lifts are sold as a kits. circled and pointed parts shows the kits that are available. Detailed list with included parts is shown at the next page.

Kit number in exploded view	Part no.	Spare parts kit	Included in kit	Quantity of pieces in kit
1a	71464	Bottom cover HeliQ - kit	Screw PT WN1452 K30x40 FZB	2
			Bottom cover	1
1b	92971	Top cover HeliQ	Cover 54028	1
			Label at customer specifications (picture of original label required)	1
2	71465	Cover plug HeliQ - kit	Cover Plug 9.3	20
3	71463	Upper limit switch Altair HeliQ - kit	Belt guide Altair HeliQ	1
			Belt stop switch	2
			Belt idler shaft	1
			Screw MC6S M6x16 A2	2
4	94356	Battery ceiling lift - kit	Battery 12V, 2.9Ah	2
			3M VHB 5930 Double sided adhesive tape	2
5	71466	Beltstop sensor Altair HeliQ - kit	Belt stop sensor with 350 cable	2
			Belt tension guard sensor 1	1
			Spring Belt Tension Guard	1
			Washer Plastic BTG	1
			Washer-LB-A-4816	1
			Belt tension guard shaft	1
			Screw MC6S M3x8 A2-70	3
			Screw MC6S M6x16 A2	2
Belt glider 1 pc.	1			
6	71477	Circuit board HeliQ 150kg - kit	Printed circuit board main HeliQ 150kg	1
			Screw PT WN1451 K30x8 FZB	2
			Washer Plastic 4.83x3.05x2.54	2
6	71478	Circuit board HeliQ 220kg - kit	Printed circuit board main HeliQ 220kg	1
			Screw PT WN1451 K30x8 FZB	2
			Washer Plastic 4.83x3.05x2.54	2
6	71479	Circuit board HeliQ 300kg - kit	Printed circuit board main HeliQ 300kg	1
			Screw PT WN1451 K30x8 FZB	2
			Washer Plastic 4.83x3.05x2.54	2
not in view	54437	Coin battery CR-1025_BN	Battery coin for main circuit board	1

Kit number in exploded view	Part no.	Spare parts kit	Included in kit	Quantity of pieces in kit
7	71490	Emergency stop HeliQ - kit	Emergency Stop HeliQ	1
8a	71483	Lift belt ceiling lift SWL 150kg/330lbs and 220kg/440lbs	Lift belt for Altair, Roomer S and HeliQ with lifting length 2120mm	1
8b	71494	Lift belt ceiling lift SWL 300kg/660lbs - kit	Lift belt for ceiling lift Altair, Roomer S, HeliQ with SWL, lifting length 1720mm	1
9	71467	Belt roll Altair HeliQ - kit	Screw MC6S M6x16 A2	6
			Belt Roll Altair HeliQ	1
			Shaft Belt Roll Altair HeliQ	1
			Washer roll	1
10a	71455	Spring telescopic hanger bar ceiling lift - kit	Spring	10
10b	72803	Locking clips for hanger bar	Clips for hanger bars 72754, 72755, 72757, 55751, 55764, 92938	10
11	71458	Hand control HeliQ - kit	Hand control 4-button HeliQ without covers	1
			Strain Relief Clamp cc12	1
			Screw PT WN1452 K30x10 FZB	2
12	71459	Hand control HeliQ Power traverse - kit	Hand control 7-button HeliQ power traverse	1
			Strain Relief Clamp cc12	1
			Screw PT WN1452 K30x10 FZB	2
13	71462	Screw holder Altair HeliQ - kit	Screw Holder	10
14	71468	Top plate screw HeliQ - kit	Screw MC6S M6x16 A2-70	20
15	71469	Pin and locking pin HeliQ - kit	Locking needle	2
			Split pin CBH 20x25 fzb	2
not in view	71470	Charger Europe, overhead lift - kit	Battery Charger	1
			Cable Charger Europe	1
			Hand control holder	1
			Ring strap for holder	1
not in view	71471	Charger Great Britain, overhead lift - kit	Battery Charger	1
			Cable Charger Great Britain	1
			Hand control holder	1
			Ring strap for holder	1
not in view	71472	Charger Japan, overhead lift - kit	Battery Charger	1
			Cable Charger Japan	1
			Hand control holder	1
			Ring strap for holder	1
not in view	71473	Charger North America, overhead lift - kit	Battery Charger	1
			Cable Charger North America	1
			Hand control holder	1

Kit number in exploded view	Part no.	Spare parts kit	Included in kit	Quantity of pieces in kit
			Ring strap for holder	1
not in view	71491	Charger Australia, overhead lift - kit	Battery Charger	1
			Cable Charger Australia	1
			Hand control holder	1
			Ring strap for holder	1
not in view	90407	Improvement color RAL9016		1
not in view	71811	Cover screw kit	Cover screw	10
not in view	92972	Secondary control panel RoomerS, Altair, HeliQ	Membrane panel Roomer S Altair HeliQ	1

Opening the lift - HeliQ

ENG

1. Activate the emergency stop.

2. Remove the top plate screws using a 5 mm Allen key. Remove the top plate.

3. Remove the two cover plugs.


4. Untighten the two screws with a torx screwdriver (TX10) holding the covers together.

5. Turn the lift around and remove the cover by angle it a bit towards the opposite side of the control panel until you hear a "click".

6. Always disconnect the battery cable when the cover is opened.

1. Make sure that all cables are connected as they were when opening the lift.
2. Make sure no cables are pinched between the covers.
3. Follow the steps 1-6 in the “opening the lift chapter” in reverse.


Don't forget the cover plugs!

Remember to attach all cables before closing the lift!

Reseting service alarm – HeliQ

After the service is finished the service alarm should be reset.

- Connect the lift to a computer thru “Data and service application”. In the software, click service and log in with your login details. Click on the button “Reset service alarm”.


The screenshot shows the 'HUMAN CARE Data- och serviceapplikation' interface. The top navigation bar includes 'START', 'ÖVERSIKT', 'STATISTIK', and 'SERVICE'. The 'SERVICE' tab is active. Below the navigation bar, there are three main sections: 'ÅTERSTÄLL SERVICEALARM' (highlighted with a red circle), 'UPPDATERA FIRMWARE', and 'STÄLL IN REALTIDSKLOCKA'. The 'SERVICE' section displays 'Senaste service' (2015-04-08) and 'Nästa serviceöppelse' (2016-04-07). The 'FELRAPPORTER' section contains a table with columns for 'Datum', 'Tid', and 'Beskrivning'.

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
The belt can be changed without removing the belt outlet.

Follow the steps 1-5 in the “opening the lift chapter”.


1. Re-activate the emergency stop.
2. Unscrew the Opto-switch for the belt tension guard and turn it 90° in order to be able to freely operate the lift without weight.


3. Drag the belt out all the way from the belt stop.

4. Create a big loop.

5. Use for example a spatula to create a loop on the other side.


6. Press “up” button on the hand control to turn the belt roll and grab the loop.

7. Use the spatula again to get the end of the belt with the small pin out.


8. Return the belt roll with help of the hand control (up-button). Grab the belt.

9. Remove the pin and the belt is loose.

9. To insert the new belt, do the steps in reverse.


After the new belt is assembled. Make sure to position the opto-switch to it's starting position!


The battery pack consists of two 12 V, 2.9 Ah, lead type batteries. The lifetime is up to 5 years depending on charging and usage.

Replacing the battery:


1. Follow the steps 1-6 in the “opening the lift chapter”.
2. Try to not damage the warning label on top of the battery, for example by folding up the ends of the label.


3. Unscrew the two metal brackets from the chassis by using a 2,5 mm Allen key.

4. Lift up the battery pack from the lift.


4. When the batteries is removed from the lift, disconnect all four cables.



5. Fit the new batteries and re-connect the cables. Make sure the batteries are connected the same way as before the change.
6. Secure the battery cables between the battery (outside the steel brackets) and the lift cover.
7. Follow the steps 1-3 in reverse to complete the battery change.


1. Follow the steps 1-6 in the “opening the lift chapter”.
2. Remove the connector on the hand control cable from the circuit board (PCB) by pressing on the small plastic lock on the left side.


3. Remove the two screws from the cable holder with a torx screwdriver 10 mm (T10).

4. Remove the rubber seal from the cover.

5. Pull the cable and the connector out through the hole.

6. Repeat the steps in reverse in order to connect the new hand control.

If the hand control is damaged/broken the cable can be cut for easier removing of the hand control and cable.

Kit number (exploded view)	Part no.	Spare parts kit	Included in kit	Quantity of pieces in kit
1a	71460	Top cover Roomer S - kit	Screw PT WN1452 K30x40 fzb	2
			Top cover Roomer S	1
1b	92936	Bottom Cover Roomer S kit	Cover 54055	1
			Label at customer specifications (picture of original label required)	1
2	71445	Foot Portable lift - kit	Foot portable lift	20
3	71446	Slider Roomer S - kit	Screw MC6LS M6x16	3
			Screw MC6S M6x16 A2	4
			Slider with sensor	2
4	71447	Switch Roomer S - kit	Roomer S switch lever	2
			Shaft Slider	1
			Screw MC6S M3x12 A2	4
			10015-N08 Mekanex	2
			Slide bearing 8x9x3	2
5	71448	Upper limit switch Roomer S - kit	Belt guide Roomer S	1
			Belt stop switch Roomer S	2
			Belt idler shaft Roomer S	2
			Screw MC6S M6x16 A2	4
6	94356	Battery ceiling lift – kit	Battery 12V, 2.9Ah	2
			3M VHB 5930 Double sided adhesive tape	2
7a	71451	Beltstop sensor Roomer S - kit	Belt stop sensor with 350 cable	1
			Belt stop sensor with 400 cable	1
			Belt tension guard sensor 2	1
			Belt tension guard sensor 1	1
			Spring belt tension guard	4
			Washer Plastic BTG	4
			Washer-LB-A-4816	2
			Belt tension guard shaft Belt tension guard shaft	2
			Screw MC6S M3x8 A2-70	4
			Screw MC6S M6x16 A2	4
7b	71551	Band tension guard Roomer S – kit, for modified chassis (after date 2021.10.05, see product label at lift)	Belt glider	2
			Cable ties 2,5x100	3
			As 7a but for ceiling lift produced after October 2021.	

Kit number in exploded view	Part no.	Spare parts kit	Included in kit	Quantity of pieces in kit
8	71452	Circuit board Roomer S - kit	Printed circuit board main Roomer S 220 kg	1
			Screw MC6S M3x8 A2-70	3
not in view	54437	Coin battery CR-1025_BN	Battery coin for main circuit board	1
9	71489	Emergency stop portable lifts - kit	Emergency Stop 16	1
10a	71483	Lift belt ceiling lift SWL 150kg/330lbs and 220kg/440lbs	Lift belt for Altair, Roomer S and HeliQ with lifting length 2120mm	1
10b	71494	Lift belt ceiling lift SWL 300kg/660lbs - kit	Lift belt for ceiling lift Altair, Roomer S, HeliQ with SWL, lifting length 1720mm	1
11a	71455	Spring telescopic hanger bar ceiling lift - kit	Spring	10
11b	72803	Locking clips for hanger bar	Clips for hanger bars 72754, 72755, 72757, 55751, 55764, 92938	10
11c	92764	Shaft hanger bar Altair, Roomer S	For assembly of hanger bar 55751 and 55764	1
12	71453	Belt roll clock wise Roomer S -kit	Screw MC6S M6x16 A2	4
			Belt roll clock wise	1
			Bushing for belt roll shaft	2
			Shaft Belt Roll	1
13	71454	Belt roll counter clock wise Roomer S - kit	Screw MC6S M6x16 A2	4
			Belt roll counter clock wise	1
			Bushing for belt roll shaft	2
			Shaft Belt Roll	1
14a	71456	Spring stainless steel leg hook Altair Roomer S - kit	Spring	5
14b	92974	Spring coated leg hook - 5 pack	Spring for leg hook	5
15	71457	Hand control Altair Roomer S - kit	Hand control portable lifts without cover	1
			Strain Relief Clamp cc12	1
			Screw PT WN1452 K30x10 FZB	2
not in view	71470	Charger Europe, overhead lift - kit	Battery Charger	1
			Cable Charger Europe	1
			Hand control holder	1
			Ring strap for holder	1
not in view	71471	Charger Great Britain, overhead lift - kit	Battery Charger	1
			Cable Charger Great Britain	1
			Hand control holder	1

Kit number in exploded view	Part no.	Spare parts kit	Included in kit	Quantity of pieces in kit
			Ring strap for holder	1
not in view	71472	Charger Japan, overhead lift - kit	Battery Charger	1
			Cable Charger Japan	1
			Hand control holder	1
			Ring strap for holder	1
not in view	71473	Charger North America, overhead lift - kit	Battery Charger	1
			Cable Charger North America	1
			Hand control holder	1
			Ring strap for holder	1
not in view	71491	Charger Australia, overhead lift - kit	Battery Charger	1
			Cable Charger Australia	1
			Hand control holder	1
			Ring strap for holder	1
not in view	90407	Improvement color RAL9016		1
not in view	71811	Cover screw kit	Cover screw	10
not in view	92972	Secondary control panel RoomerS, Altair, HeliQ	Membrane panel Roomer S Altair HeliQ	
not in view	92973	Battery holder kit	54905 Battery holder RoomerS	1
			55629 Heat shrink tubing	1

Opening the lift - Roomer S

ENG

1. Activate the emergency stop.


2. Turn the lift upside down. Be careful with the belt stop sensor.

3. Remove the rubber feet closest to the hanger bar on Roomer S lift.


4. Untighten the two screws (TX10) holding the covers together. The screws are located at the bottom of the holes for the rubber feet.

5. Turn the lift around and remove the cover by angle it a bit towards the opposite side of the control panel until you hear a "click".

6. Always disconnect the batteries when the cover is opened.

1. Make sure that all cables are connected as they were before.
2. Make sure no cables are pinched.
3. Follow the steps 1-6 in the “opening the lift chapter” in reverse.


Remember to attach the battery cable before closing the lift!

Don't forget the rubber feet!

Resetting service alarm – Roomer S

After the service is finished the service alarm should be reset.

- Connect the lift to a computer thru “Data and service application”. In the software, click service and log in with your login details. Click on the button “Reset service alarm”.

The screenshot shows the 'Data- och serviceapplikation' interface for HUMAN CARE. The top navigation bar includes 'START', 'ÖVERSIKT', 'STATISTIK', and 'SERVICE'. The 'SERVICE' tab is active. Below the navigation bar, there are three main sections: 'ÅTERSTÄLL SERVICEALARM' (circled in red), 'LIPPOÄTTRA TIDPUNKTER', and 'STÄLL IN REALTIDSKLOCKA'. The 'SERVICE' section displays a table with columns for 'Senaste service' and 'Nästa service tillfälle'. The 'FELRAPPORTER' section displays a table with columns for 'Datum', 'Tid', and 'Beskrivning'.

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


Changing the hand control - Roomer S

ENG

1. Follow the steps 1-6 in the “opening the lift chapter”.
2. Remove the connector on the hand control cable from the circuit board (PCB) by pressing on the small plastic lock on the top.


3. Remove the rubber seal from the cover.


4. Use for example a screwdriver to flip over

the cable to the left side like in pictures below.


5. Pull the cable and the connector out through the hole. If the hand control is damaged, you can cut the cable and it will be easier to pull out.


6. Repeat the steps in reverse in order to connect the new hand control.

The belt can be changed without removing the belt outlet.


1. Drag the unlocked belt out all the way from the belt roll counter.


2. Create a quite big loop.

3. Use for example a spatula to create a loop on the other side.

5. Use the roomer S switch or the hand control to turn the belt roll and grab the loop.

6. Use the spatula again to get the end of the belt with the small pin out.


7. Return the belt roll with help of the Rooter S switch or the hand control. Grab the belt.

Note! Do not turn the Rooter Switch without holding the belt roll since the spring will unwind. If you drop the belt roll, rewind the belt roll 9,5 turns from the belt roll slot with a metal pin.


7. Remove the pin and the belt is loose.

8. If you need to remove both belts, do the same on the opposite side.

9. Insert the new belt/belts through the belt outlet as shown in the picture below.

Note! Make sure the black tags on both belts are facing each other.

Note! Make sure the belts are on the correct side of the Belt tension guard.


- The belt should enter the belt roll through the slot with the metal pin.
- The end of the belt should be fastened with the metal pin in the slot next to the slot where it entered the belt roll.
- The pin should be resting in the slot as much as possible. Pull the belt to make the pin fit in the slot.

The battery pack consists of two 12 V, 2.9 Ah, lead type batteries. The lifetime is up to 5 years depending on charging and usage.

Replacing the battery:


1. Follow the steps 1-6 in the “opening the lift chapter”.

Note! Make sure to remove the battery cable that connects to the circuit board (not the battery cables that is placed on the battery).


2. Unscrew the metal bracket from the chassis by using a 2,5 mm Allen key.

Warning! Do not remove the battery cables show in the picture below! They should stay on until the screw is loosened to avoid short circuit.


3. Now - disconnect the battery cables located on the battery.

4. Lift up the battery pack from the lift.


5. When the batteries is removed from the lift, disconnect all four cables.

5. Fit the new batteries and re-connect the cables. Make sure the batteries are connected the same way as before. The batteries connectors from the circuit board and the battery terminals are color coded. Cables from emergency stop can be connected in any order.
6. Secure the battery cables between the battery (outside the steel brackets) and the lift cover.
7. Follow the steps 1-5 in reverse to complete the battery change.

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Kitnummer i sprängskiss	Art. nr.	Reservdelskit	Inkluderat i kit	Antal delar i kit
1a	71461	Övre kåpa Altair - kit	Screw PT WN1452 K30x40 FZB	2
			Top cover Altair	1
1b	92937	Undre kåpa Altair - Kit	Cover 54027 Label at customer specification (picture of original label required)	1 1
2	71445	Fötter portabel lyft - kit	Foot portable lift	20
3	71463	Övre brytare Altair HeliQ - kit	Belt guide Altair HeliQ	1
			Belt stop switch	2
			Belt idler shaft	1
			Screw MC6S M6x16 A2	2
4	94356	Batteri taklyft - kit 2 st	Battery 12V, 2.9Ah	2
			3M VHB 5930 Double sided adhesive tape	2
5	71466	Sensorer bandstyrning Altair HeliQ - kit	Beltstop sensor with 350 cable	2
			Belttension guard sensor 1	1
			Spring Belt Tension Guard	1
			Washer Plastic BTG	1
			Washer-LB-A-4816	1
			Belt tension guard shaft	1
			Screw MC6S M3x8 A2-70	3
			Screw MC6S M6x16 A2	2
6	71474	Kretskort Altair 150kg - kit	Printed circuit board main Altair 150kg	1
			Screw PT WN1451 K30x8 FZB	2
6	71475	Kretskort Altair 220kg - kit	Printed circuit board main Altair 220kg	1
			Screw PT WN1451 K30x8 FZB	2
6	71476	Kretskort Altair 300kg - kit	Printed circuit board main Altair 300kg	1
			Screw PT WN1451 K30x8 FZB	2
7	71489	Nödstopp portabel lyft - kit	Emergency Stop 16	1
8a	71483	Lyftband taklyft SWL 150kg och 220 kg- kit	Lyftband för Altair, Roomer S, HeliQ, med lyftlängd 2120 mm	1
8b	71494	Lyftband taklyft SWL 300 kg - kit	Lyftband för Altair, Roomer S, HeliQ, med lyftlängd 1720 mm	1
9	71467	Bandrulle Altair HeliQ - kit	Screw MC6S M6x16 A2	6
			Belt Roll Altair HeliQ	1
			Shaft Belt Roll Altair HeliQ	1
			Washer roll	1

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Kitnummer i sprängskiss	Art. nr.	Reservdelskit	Inkluderat i kit	Antal delar i kit
1a	71460	Övre kåpa Roomer S - kit	Screw PT WN1452 K30x40 fzb	2
			Top cover Roomer S	1
1b	92936	Undre kåpa Roomer S kit	Kåpa 54055	1
			Etikett vid kundspecifikation (bild av original etikett nödvändig)	1
2	71445	Fötter portabel lyft - kit	Foot portable lift	20
3	71446	Glidplatta Roomer S - kit	Screw MC6LS M6x16	3
			Screw MC6S M6x16 A2	4
			Slider with sensor	2
4	71447	Spak Roomer S - kit	Roomer S switch lever	2
			Shaft Slider	1
			Screw MC6S M3x12 A2	4
			10015-N08 Mekanex	2
			Slide bearing 8x9x3	2
5	71448	Övre brytare Roomer S - kit	Belt guide Roomer S	1
			Belt stop switch Roomer S	2
			Belt idler shaft Roomer S	2
			Screw MC6S M6x16 A2	4
6	94356	Batteri taklyft - kit 2 st	Battery 12V, 2.9Ah	2
			3M VHB 5930 Double sided adhesive tape	2
7	71451	Sensorer bandstyrning Roomer S - kit	Beltstop sensor with 350 cable	1
			Beltstop sensor with 400 cable	1
			Belttension guard sensor 2	1
			Belttension guard sensor 1	1
			Spring belt tension guard	4
			Washer Plastic BTG	4
			Washer-LB-A-4816	2
			Belt tension guard shaft	2
			Screw MC6S M3x8 A2-70	4
			Screw MC6S M6x16 A2	4
			Cable ties 2,5x100	3

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The screenshot shows the 'HUMAN CARE' Data- och serviceapplikation interface. The 'SERVICE' tab is selected. The 'ÅTERSTÄLL SERVICEALARM' button is circled in red. Below it, there are buttons for 'UPPDATERA FIRMWARE' and 'STÄLL IN REALTIDSKODEN'. The 'SERVICE' section shows the latest service date as 2015-04-08 and the next service due date as 2016-04-07. The 'FELRAPPORTER' section contains a table of error reports.

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