



Human Care's new freestanding gantry FlexiRail is a perfect alternative to ceiling-mounted rails. FlexiRail works perfectly with all our lifts and can be adjusted laterally and in height. It can also be equipped with wheels for easy movement of the system between rooms.

#### Portable Solution

Sometimes a facility or homeowner does not want to make changes to the ceiling by installing a permanent rail system. For example, in the event of a temporary need for a lifting solution. In such a case our independent rail system FlexiRail is the perfect choice. Portable, easy to use and to assemble without the need for interventions to walls or ceilings. FlexiRail's compatibility with our ceiling lifts Altair, HeliQ and Cosmos makes it possible to use FlexiRail for many patients and in many different lifting situations.

#### Flexibility at its best

FlexiRail is available in several different lengths, and as it is adjustable with a range between 265-475cm / 104.3-187", FlexiRail

can be completely adapted to your needs. The stand can also be adjusted in height, from 180-270cm / 70.9-106.3". The adjustment of height and length takes place easily without the need for tools. In combination with our rail rollers, ceiling lifts, hanger bars and slings, a customized lifting solution can be created. FlexiRail blends aesthetically into most environments thanks to the powder-coated white surface.

#### Quality

For us, safety, quality and usability are of the utmost importance. We constantly strive to ensure and improve the highest quality of all our products. In every detail, in every aspect, our products keep what they promise.



## Accessories



### Wheels

Optional wheel kit to assist with portability. The height of the wheel including the mounting plate is 10 cm/3,9". Part no. 50872

**Human Care offers a wide range of rail rollers, ceiling lifts, hanger bars and slings as well as other lifting accessories for a customized solution.**

## Lifting Situations



### Use to/from a bed

Can be placed over bed for transfers and repositioning with a ceiling lift/sling.



### Use to/from a chair

Can be placed over a chair for transfers with a ceiling lift/sling.



### Use to/from the floor

Can be placed over a patient in the supine position for transfers/lifting with a ceiling lift/sling.

## Technical specifications

Product name/ Article number	Height cm / "	Length cm / "	Width, feet cm / "	SWL kg / lb
FlexiRail (250 cm) / 98.4" 90900-250	180 - 270 cm 70.9 - 106.3"	265 - 325cm 104.3 - 127.9"	120 cm 47.2"	300kg / 660lb
FlexiRail (300 cm) / 118.1" 90900-300	180 - 270 cm 70.9 - 106.3"	315-375cm 124 - 147.6"	120 cm 47.2"	300kg / 660lb
FlexiRail (350 cm) / 137.8" 90900-350	180 - 270 cm 70.9 - 106.3"	365-425cm 143.7 - 167.3"	120 cm 47.2"	300kg / 660lb
FlexiRail (400 cm) / 157.5" 90900-400	180 - 270 cm 70.9 - 106.3"	415-475cm 163 - 187"	120 cm 47.2"	300kg / 660lb

98054\_Rev1\_EN\_Product sheet\_FlexiRail\_Aug24



*We improve freedom of movement  
and quality of care for people.*